

AccuTOF-GCv Series

Analysis of a Block Copolymer by Using Field Desorption (FD)

Introduction

A commercially available polyoxypropylene (PO) - polyoxyethylene (EO) block copolymer was analyzed by using the JMS-T100GC "AccuTOF GC" field desorption (FD) method. A group-type analysis was performed on the resulting mass spectrum.

Methods

Sample

Poly(ethylene glycol)-*block*-poly(propylene glycol)-*block*-poly(ethylene glycol) (Aldrich 435406; BASF Pluronic® L-31), $M_n = 1100$, 50 mg/mL in THF




Fig. 1 Structural formula of the sample

MS conditions

Mass spectrometer: JMS-T100GC "AccuTOF GC"
 Ionization mode: FD (+)
 Cathode potential: -10 kV
 Emitter current: 0 mA \rightarrow 51.2 mA/min \rightarrow 40 mA
 Acquired mass range: m/z 35 – 1,600
 Spectral recording interval: 1.0 sec

Results and discussion

The FD mass spectrum for the block copolymer sample is shown in Fig. 2. The base peak was observed at m/z 1005.7 with the other dominant peaks (m/z 947.7, m/z 889.7, etc.) separated by an interval of 58Da between them that corresponds to the loss of $\text{C}_3\text{H}_6\text{O}$. These ions represent the structure shown in Fig. 1 where $x = 0$, $y = 15$, 16, 17, and $z = 0$.

For group-type analysis, the following parameters were used:

Group-type analysis parameters

Software: Polymerix™ (Sierra Analytics)
 Repeat unit A: $\text{C}_3\text{H}_6\text{O}$
 Repeat unit B: $\text{C}_2\text{H}_4\text{O}$
 α end group: H
 ω end group: OH
 Adduct: H
 Match tolerance: ± 0.05 u




Fig. 2 FD mass spectrum





Fig.3 Type analysis assignments (m/z 920-1020)



Despite the fact that the sample was an EO-PO-EO triblock copolymer, only two repeat units, PO and EO, were specified as the repeat units since two of the three blocks were EO. The peak assignments for the m/z 920-1020 range are shown in Fig. 3 where “A” is C_3H_6O , “B” is C_2H_4O , and the number following “A” or “B” represents the number of repeat units. The peaks that correspond to $A = 12$ to 17 and $B = 0$ to 10 have been assigned in Fig. 3. The distribution of the A (PO) and B (EO) repeat units are shown in Fig. 4. Also, the average molecular weights and other metrics derived from the group-type analysis are shown in Table 1.

Conclusions

By analyzing the FD mass spectrum of a triblock copolymer with a suitable group-type analysis software, metrics such as ratios and distribution among constituents, number average molecular weight (M_n), weight average molecular weight (M_w), Z average molecular weight (M_z), and polydispersity (PD) were easily obtained using the AccuTOF GC field desorption method.

		Copolymer Distribution												
		Repeat B												
		0	1	2	3	4	5	6	7	8	9	10	11	
Repeat A	6													
	7													
	8													
	9						1.08	1.66	1.27	1.33	1.61			
	10			1.63	2.17	4.19	4.24	3.04	2.26	2.58	1.72			
	11	2.65	3.42	5.09	6.17	8.20	7.60	5.81	4.63	4.34	2.42			
	12	9.58	10.37	12.54	14.00	15.99	13.09	9.23	5.62	4.32	3.27	2.19	1.40	
	13	25.98	23.80	23.46	24.20	23.93	17.79	12.39	7.36	6.81	3.56	2.24	1.22	
	14	54.36	39.41	35.20	31.16	26.41	20.61	13.98	9.37	7.03	4.81	2.82	1.11	
	15	83.18	52.75	43.08	34.49	28.19	19.84	14.31	9.97	6.35	2.97	2.14	1.09	
	16	99.99	56.93	41.44	32.28	25.35	18.00	12.60	8.69	6.50	3.50			
	17	100.00	52.72	36.53	26.54	20.38	14.20	10.73	6.67	3.61				
	18	84.89	41.17	27.92	20.53	16.96	9.92	7.72	4.55	3.09	1.37			
	19	63.64	29.20	19.25	14.66	10.89	7.26	4.05	3.16	1.73				
	20	40.33	19.86	13.74	9.46	6.85	4.96	1.56	1.69					
	21	25.17	11.41	8.93	5.77	4.50	3.04	2.11						
	22	14.63	7.35	3.67	3.55	2.94								
23	7.94	4.08	1.86											
24	4.44	2.56												
25	1.69													
26														

Fig.4 Copolymer distribution

	M_n	M_w	M_z	PD
$(H[C_3H_6O]_n[C_2H_4O]_mOH) + H^+$	1053.2	1077.1	1100.9	1.03

Table 1. Type analysis result