

JEOL

MULTIBeam


MultiBeam

High performance SEM and micro milling FIB

New, easy-to-use MultiBeam system combines high resolution ion optics with the most popular SEM column in the world. The JIB-4500 MultiBeam offers increased throughput and productivity for a variety of applications, from viewing to analysis to micro milling.


Large current milling mode

- Enhances throughput from large milling area (maximum milling current 30 nA).



High resolution electron (LaB₆) and ion optics

- Provides stable imaging for milling/monitoring over extended periods of time.



Electron column resolution


Ion column resolution

Charge-free imaging and analysis

- Low vacuum system standard for non-conductive specimens

Bitmap milling system

- Enables deposition and milling of binary image data stored in BMP file.
- Supports 3D nano milling with continuous access to multiple bitmap data.



Standard airlock system for fast and easy sample loading

- Fast sample loading through airlock chamber.
- Maximum sample size: 150mm wafer.

Low mag mode identifies milling positions

- Milling position initially identified at minimum magnification of 5X in SEM.
- Transition to milling mode is easily made at minimum magnification of 30X in FIB.

MultiBeam *High performance SEM and micro milling FIB*




Monitor milling in real time

- Simultaneous SEM imaging of FIB milling process.
- Effectively prepares TEM thin film samples to optimum thickness.

A complete set of ports for all your applications


- Large stage is capable of image milling of sample areas up to 76mm.
- A complete set of ports for a variety of analyses, including EDS, EBSD, and CLD, are arranged at optimum positions for streamlining workflow.



* Optional accessories


Optional chamberscope verifies sample positions

- Effective for verifying positional relationships with other options including gas injection system.




Serial slicing and sampling - S³™

- Monitor slicing and fabrication in process.
- Reconstruct 3D images



User friendly, easy to learn

- Modern, clearly marked GUI.
- Visual, intuitive operation from SEM/FIB condition setting to milling and monitoring.



MultiBeam

High performance SEM and micro milling FIB

Principal Specifications

FIB

Ion source	Ga liquid metal ion source
Resolution	5nm @ 30kV
Accelerating voltage	1 to 30kV (in steps)
Magnification	30x (for wide view) 100X to 300,000X
Maximum beam current	30nA (at 30kV)


SEM

Beam source	LaB6
Resolution	2.5nm @ 30kV
Accelerating voltage	0.3 to 30kV (in steps)
Magnification	5X to 300,000X
Maximum beam current	1,000nA (at 30kV)
Goniometer stage	X: 76mm; Y: 76mm; Z: 5 to 48mm Z' for eucentric point adjustment T: -10° to 90°; R: 360° endless
Vacuum pump	SIP (x2), TMP/RP (x2)
Standard features	Low vacuum (LV) Backscatter detector (BEI) Gas injection unit (1)

Installation Requirements


Power	Single phase 100V± 10% 50/60 Hz, 3 kVA (x2)
Ground	One, 100 ohm or less
Dry nitrogen gas	0.4 to 0.7 MPa
Room temperature	18 to 25°C
Humidity	60% or less

Specifications subject to change without notice.



Optional Attachments

Chamber scope - top mount/side mount
Gas injection system
External scan interface for SEM
EDS, EBSD, CLD



JEOL USA, Inc.

11 Dearborn Road • Peabody, MA 01960
Tel: 978-535-5900 • FAX: 978-536-2205

For more information, contact your local sales office,
email salesinfo@jeol.com, or visit www.jeolusa.com.

5653 Stoneridge Drive • Pleasanton, CA 94588
Tel: 925-737-1748 • FAX: 925-737-1794