

Elemental Compositions from Exact Mass Measurements and Accurate Isotopic Abundances

Introduction

Exact masses have been used for decades to calculate elemental compositions for known and unknown molecules. The traditional approach calculates all possible combinations of user-specified atoms that fall within a given error tolerance of a measured mass. The number of possible combinations increases dramatically with increasing mass and as more atoms are included in the search set. In many cases, it is not possible to determine a unique composition based on mass alone¹.

A common source of error in measuring isotopic abundances with scanning mass spectrometers is

related to fluctuations in ion current during measurement. The AccuTOF family of mass spectrometers overcomes this problem by analyzing all of the isotopes formed at the same instant. Combined with a high-dynamic-range detector, this provides highly accurate isotopic abundances². It has been shown that accurately measured isotopic abundances can be combined with measured exact masses to dramatically reduce the number of possible elemental compositions for an unknown³. It is often possible to deduce a unique elemental composition, facilitating the identification of unknown substances.




Figure 1. Elemental composition calculation for thioridazine from combined exact mass measurement with isotope matching.

Experimental

Samples in this report were measured with the AccuTOF-DART[™] mass spectrometer. Similar procedures can be used with other members of the AccuTOF mass spectrometer family. Calibrated mass spectra were centroided and saved as JEOL-DX (JCAMP) text files. These text files were processed with the Elemental Composition Workshop from the Mass Spec Tools[™] software suite distributed with AccuTOF mass spectrometers⁴. The program permits automated isotope matching for measured mass spectra and provides a visual comparison between the measured and theoretical isotopic abundances for each hit (Figure 1).

Eleven small drug samples were deposited on melting point tubes and mounted on a support. Samples were passed sequentially in front of the DART source with a measurement time of about 3 seconds per sample. Neat PEG 600 on a melting point tube was also measured to provide an external

mass calibration standard. Elemental compositions were calculated by assuming even-electron ions, a mass measurement error tolerance of 0.003 u, and default tolerances for isotope matching. Elemental limits were set to:

C 0/70 H 0/80 O 0/20 N 0/10 S 0/2 Cl 0/1.

The correct elemental composition was successfully determined for each sample as shown in Table 1. The total number of compositions calculated without isotope matching is shown in column 4 of the table. The correct composition was the number 1 ranked composition from the automated isotope match for all compounds measured. The most dramatic example was reserpine [M+H]⁺ (C₃₃H₄₁N₂O₉) which gave 30 compositions without isotope matching. Isotope matching reduced the number of possible compositions to 11 of which the best match corresponded to the correct composition.

Compound	Composition	Calculated <i>m/z</i>	#Compositions	Rank
Phenolphthalein	C ₂₀ H ₁₅ O ₄	319.097035	18	1
Promazine	C ₁₇ H ₂₁ N ₂ S	285.142544	11	1
Quinine	C ₂₀ H ₂₅ N ₂ O ₂	325.191603	8	1
Nortriptylene	C ₁₉ H ₂₂ N	264.175224	5	1
Thioridazine	C ₂₁ H ₂₇ N ₂ S ₂	371.161565	16	1
Chlorpromazine	C ₁₇ H ₂₀ N ₂ SCl	319.103572	18	1
Doxepin	C ₁₉ H ₂₂ NO	280.170139	5	1
Hydrocortisone	C ₂₁ H ₃₁ O ₅	363.21715	9	1
Reserpine	C ₃₃ H ₄₁ N ₂ O ₉	609.281208	30	1
Caffeine	C ₈ H ₁₁ N ₄ O ₂	195.088201	7	1
Erythromycin	C ₃₇ H ₆₈ NO ₁₃	734.469069	22	1

Table I. [M+H]⁺ elemental compositions determined for 11 drugs measured in 0.58 minutes with the AccuTOF-DART. All compositions were correctly identified. The rms mass measurement error was 2.2 ppm.

References

1. <http://jeolusa.com/ms/docs/elcomp.html>
2. <http://jeolusa.com/ms/docs/phenylalanine.pdf>
3. <http://www.epa.gov/nerlesd1/chemistry/ice/pubs.htm>
4. JEOL USA, Inc. customers only.

The software is available to all for purchase at <http://www.chemsw.com/13057.htm>