

AccuTOF-GCv Series

Group-type Analysis of Crude Oil by Using GC/FI-TOFMS 2

Reproducibility of group-type analysis results

Introduction

Field Desorption (FD) and Field Ionization (FI) are both techniques that ionize analytes by electron tunneling from the analyte molecules to a solid surface (emitter) in a high electric field. In the case of FD, the sample is applied directly onto the emitter and heated by applying an electric current through the emitter for desorption and ionization. In the case of FI, vaporized analyte molecules are introduced into the proximity of the emitter.

Both FI and FD are soft ionization methods that generally yield intact molecular ions and, in most cases, produce very few fragment ions. Generally, these two techniques are used to ionize analytes that are easy to fragment and do not generally produce molecular ions during electron ionization (EI), such as hydrocarbons in crude oil.

For complex mixtures such as crude oils or synthetic polymers, molar mass distributions and average molecular weights are important chemical properties. By analyzing an FI mass spectrum that mostly consists of molecular ions, even for complex hydrocarbon mixtures, and using a group-type analysis software, one can obtain the molar distributions and average molecular weights for the various hydrocarbon types (e.g., paraffin, naphthene, olefin, aromatics) present in the sample mixture.

In this work, we analyzed a crude oil sample 10 times by using the JMS-T100GC "AccuTOF-GC" GC/FI method and then processed the resulting data using a group-type analysis software package. Afterwards, these results were then compared to each other to confirm the reproducibility of the analyses.

Method

Sample Crude Oil
(Refer to Fig. 1 for Sample Preparation)

GC Conditions

GC: Agilent 6890N
 Column: DB-5ms
 30 m x 0.25 mm I.D., 0.25 μ m
 Oven: 50°C \rightarrow 15°C/min \rightarrow 280°C
 (5 min)
 Injection port: 280 °C, Split (1:200)
 Injection volume: 1.0 μ l
 Carrier gas: He (1mL/min, constant flow mode)

MS Conditions

MS: JMS-T100GC "AccuTOF GC"
 Ionization: FI+ (Cathode voltage: -10 kV,
 Emitter current: 0 mA)
 Mass range: m/z 35 - 500
 Acquisition rate: 0.3 s/spectrum

Software

Polymerix™ (Sierra Analytics, Inc.)




Fig. 1 Sample preparation flow

