



# AccuTOF-GCx Series

## Introduction of latest GC/HR-TOFMS system: JMS-T200GC AccuTOF GCx - High Mass Resolution

### Introduction

JEOL has recently announced a 4<sup>th</sup> generation GC/HR-TOFMS system, the JMS-T200GC AccuTOF GCx, in 2015. The AccuTOF-GCx offers high sensitivity, high mass resolving power, high mass accuracy, and a wide dynamic range in combination with high-speed data acquisition.

In this application note, we show the high mass resolution that can be obtained with this latest GC/HR-TOFMS system.

### Results and Discussion

We measured perfluorotributylamine (PFTBA) introduced through the reservoir of the heated volatiles inlet. A resolving power of 10,000 was easily obtained for the  $m/z$  614 peak (Figure 2).

Next, we checked the mass separation for a doublet comprised of the peaks corresponding to  $C_6F_{11}^+$  (PFK) and  $C_7H_{21}O_4Si_4^+$  (GC column background). In comparison with the original 1<sup>st</sup>-generation model, the “AccuTOF GC” (with only a mass resolving power of 5,000), the 4<sup>th</sup>-generation model “AccuTOF GCx” has over 10,000 mass resolving power (Figure 3).



Figure 1. 4<sup>th</sup> generation of the JEOL GC/HR-TOFMS system JMS-T200GC AccuTOF GCx

### Conclusion

The AccuTOF GCx system has a mass resolving power of  $>10,000$  that can easily be obtained for all GC/MS and direct MS measurements.




Figure 2. PFTBA EI mass spectrum




Figure 3. Evolution of mass resolving power from the 1<sup>st</sup>-generation AccuTOF-GC to the 4<sup>th</sup>-generation AccuTOF-GCx