

JCM-7000, NeoScope™ Benchtop SEM


Live 3D

REAL-TIME 3D SURFACE RECONSTRUCTION WITH NEOSCOPE™

SEM is a natural extension to viewing specimens with an optical microscope due in part to its inherent higher depth of field and ability to resolve smaller microstructures. Creating a 3-dimensional (3D) surface model can further enhance our understanding with specimens that have complex topographical features.



The JCM-7000 includes a high sensitivity, multi-segmented backscatter electron detector that allows for collection of multiple images simultaneously. These images are then combined to provide a 3D model of the surface in Real-Time as you navigate around the specimen. The Live 3D surface rendering can be tilted and rotated and the resulting 3D image output directly to a report.



Live 3D Image - Pharmaceutical Tablet Surface

With the addition of JEOL's Smile View™ Map software not only qualitative but also quantitative texture information can be obtained and much more.

Use Smile View™ Map to:

- Improve image quality
- Colorize SEM images
- Apply texture analysis to the 3D surface models
- Stitch images
- Correlate with data collected from other instruments or detectors

ISO 25178	
Height Parameters	
Sq	8.14 μm
Ssk	-0.933
Sku	2.76
Sp	12.6 μm
Sv	21.5 μm
Sz	34.1 μm
Sa	6.68 μm

