

# JEM-9320FIB

*Focused Ion Beam System for Specimen Preparation*



**JEOL**

*Serving Advanced Technology*

# JEM-9320FIB

## Focused Ion Beam System for Specimen Preparation

Nanoscale observation and analysis are increasingly important in the nanotechnology era. Based on electron- and ion-optics technologies developed by JEOL for many years, the JEM-9320FIB system can prepare cross sections from local areas of STEM/TEM specimens for semiconductor failure analysis and achieves precise, fast milling for SEM cross sectional specimens.

### Preparation of SEM cross sectional specimens

- For SEM common specimen stub  
(When a bulk specimen stage is installed)  
After cross sectional milling, it is possible to transfer a specimen, which is mounted on a specimen stub, from the JEM-9320FIB to a JEOL SEM, EPMA or Auger instrument for observation and analysis.

### Simple thin-film preparation for STEM and TEM

- For TEM common tip-on holder  
(When a side-entry goniometer stage is installed)  
Observation using a JEOL TEM and re-processing using this system can be easily repeated.

### High-speed, high-precision milling, high-resolution observation

- High-speed milling with a large ion-beam current (30 kV, 30 nA or greater) reduces the rough-milling time to 1/8 (compared to JEOL conventional models).
- High-speed processing with a large ion-beam current (30 kV, 30 nA or greater) reduces the rough-processing time to 1/8 (compared to JEOL conventional models).
- Improvement of the ion optical system enables one to observe a high-resolution SIM image (6 nm, 30 kV).
- Stable beam scan allows high-precision processing.
- The optional Automatic Processing Software has greatly increased the processing efficiency.

### Eco-friendly, space-saving, energy-saving


- The JEM-9320FIB can be installed even in a narrow space because of its small footprint.
- The beam saver mode suppresses unnecessary consumption of the Ga ion source.

\* SIM image (Scanning Ion Microscope Image)

# Simple, Fast, Precise for Specimen Prep

## SEM stage/SEM common specimen stub

It is possible to use a JEOL SEM, EPMA or Auger instrument with a common specimen stub. A milled specimen (mounted on the specimen stub) by FIB can be transferred to an instrument with the specimen stub attached, for observation and analysis. Since the same specimen stub can be used, the specimen can be easily positioned for an additional processing point.




(Optional attachments are included in the photograph of the appearance of the JEM-9320FIB.)

## TEM stage/TEM common tip-on holder

The JEM-9320FIB employs a side-entry goniometer stage commonly used with a JEOL TEM. Since a tip-on holder and a shuttle retainer can be shared with a JEOL TEM, FIB milling and TEM observation can be easily repeated.

\* For the compatibility with a TEM, please consult us.




FIB tip-on holder with a shuttle retainer attached to the tip

# se. FIB System Optimum aration




SIM image

[Specimen: 512 DRAM]









SEM image

[Specimen: 512 DRAM]



STEM image/EDS maps

[Specimen: 512 DRAM]




# Twin Stage Optimum for SEM/TEM Specimen Preparation

The JEM-9320FIB can simultaneously install the bulk specimen motor stage optimum for SEM/TEM specimen preparation from a bulk specimen, and the side-entry goniometer stage optimum for TEM specimen preparation (either one is a standard configuration). A better one is selectable depending on specimen size.

## Bulk specimen stage


The entire surface of a specimen with 20 mm × 20 mm can be observed.



Bulk specimen motor stage

## Airlock system for simple specimen exchange (optional)

A specimen holder can be installed on the specimen stage through an airlock chamber in a short time.



Airlock system

## TEM specimen stage

It is possible to prepare a thin film from a specimen polished by dicing in prior treatment.




Side-entry goniometer stage for TEM specimens

## Tip-on holder for simple handling of specimens (optional)


The shuttle retainer can be installed on the tip of the tip-on holder. Since the retainer can be replaced with the specimen mounted on it, even a small specimen 3 mm in width can be easily handled.



Shuttle retainer



Shuttle retainer installed on the tip of the FIB tip-on holder



## FIB bulk specimen holder (optional)


This holder can retain a bulk specimen up to 8 mm in diameter and 1 mm thick. It is useful for cutting a thin section from a fragment of a semiconductor wafer for FIB milling.

# High Operability and Extensive Capabilities

## Wide view, snap shot


The wide view mode facilitates the search of target-milling positions by means of an ultra-low magnification ( $\times 50$ ). In particular, this capability is useful when you want to find target milling-positions on a large specimen and multiple specimens. The snap shot can capture a real-time image on the monitor screen. By clicking on an arbitrary position on a snap shot, this specified position appears at the center of the real-time screen. Combined use with the wide view mode opens the application for navigation.

Clicking the snap shot screen displays a high-magnification image. A sophisticated eucentric capability maintains the target specimen position on the monitor screen even the specimen is rotated or tilted. Sequentially varied magnifications are displayed from  $\times 150$ .




SIM image at a magnification of  $\times 150$

Wide view mode ( $\times 50$ )/Wide field image (3 mm in width)




Facilitated rotation and tilting by eucentric capability



## Thumbnail position-memory

This capability makes it possible to memorize the presently displayed image and to display it in the thumbnail view, which includes the positional information. Clicking an image in the thumbnail view moves the specimen stage and reproduces the target object at the center of the screen.


It is easy to reproduce the target specimen position when multiple specimens are alternately observed or additional milling is applied to the specimen after SEM observation. In addition, this advantage is useful when the specimen is removed from the specimen stage due to an interruption work and returned to the stage for the observation.



## Specimen position display and specimen rotation capability


The specimen position on the stage is graphically displayed. Clicking on the graphic display moves the stage to the target position and its image appears. This capability is useful for searching the target position.

With a one-touch operation on the displayed image, the specimen can be rotated to an arbitrary angle. It is easy to adjust the orientation of the specimen for processing.





## Automatic Processing Software (optional)

This software records the processing procedure (recipe) for metal deposition, rough processing and final processing in advance. In addition, the positions to be processed can be recorded up to 250 points for continuous automatic processing. This capability is utilized for all-night automatic processing and operators' final processing in the next morning, thus increasing efficiency.




## Large-current processing

Large-current processing (30 nA or more) has greatly shortened the processing time. The rough-processing time has been reduced to 1/8 (compared to JEOL conventional models). This capability is particularly useful for processing on large areas.



Cross-section SIM observation  
Specimen: bonded wire.




Large-area processing in the large-current mode  
(rough processing: 40 nA).


# Applications of JEM-9320FIB To Processing and Observation

## Cross-section processing/SIM image observation

The JEM-9320FIB is an effective tool for not only processing but also observation. In particular, SIM images offer better channeling contrast than SEM images, which are formed by differences in crystalline orientations. This advantage is very suited to metallic textures and plated films.




SIM image of a chip capacitor on a printed circuit board



SIM image of a cross section of a chip capacitor milled by an ion beam. Channeling contrast clearly shows the difference in layer structures.


## Specimen preparation for SEM

A cross section milled by FIB can be subject to high-resolution SEM observation and elemental analysis using an energy dispersive X-ray spectrometer.




SIM image

Specimen: 512M DRAM



SEM (secondary electron) image


## Elemental mapping images obtained by EDS



Composition images formed by backscattered electrons

## Thin-film specimen preparation


The JEM-9320FIB can prepare thin films from local microscopic areas with high accuracy. This tool is essential for evaluation and defect analysis of nano-materials such as semiconductor devices.




## Fine fabrication

The JEM-9320FIB can fabricate the tip of a cantilever as well as prepare thin films.

### Tip of a cantilever for an AFM fabricated by an ion beam



### Cantilever for $\alpha$ -SNOM subject to cross-section fabrication




## Thinning by additional processing


Additional processing is effective for preparation of a thin film when the area of observation is specified. In this technique, first, an ion beam cuts a specimen into a block 10 to 15  $\mu\text{m}$  in size, including the area of observation. Then, this block is picked up with a glass probe, fixed on a specimen mount such as a half-sliced grid, and attached to a shuttle retainer. This shuttle retainer is installed in the FIB system and the specimen is thinned for TEM/STEM observation. Since TEM observation and FIB processing can be repeated, the target position is precisely thinned.

### Light-emitting diode (LED) milled with a Cross Section Polisher for cross-section preparation


A LED was subject to cross-section preparation with a Cross Section Polisher for SEM observation. After this, the LED was cut into a block by FIB, fixed on a specimen mount, and thinned by FIB for high-resolution STEM observation.




Composition image formed by backscattered electrons taken with an FE-SEM



Specimen block cut by an ion beam (interface of gold bonding and LED shown in SIM image)





STEM image of a picked-up thin film



Thinning by FIB (SIM image)

It is easy to repeat thinning until an optimum specimen thickness is obtained.



STEM image of an interface and elemental mapping images obtained by EDS

# Optional Attachments for Supporting Specimen Preparation

## Specimen pickup system

This system picks up a TEM specimen in the air, which is milled by the JEM-9320FIB, and transfers it to a grid for observation. The system is comprised of a micropipette preparation device, a glass probe preparation device and an optical microscope with a manipulator.



Micropipette preparation device




Glass probe preparation device




Optical microscope with manipulator  
(CCD camera is separately provided)


## Examples of applications (thin-film pickup)



Thin film cut and separated




Thin film picked up by  
a glass probe



STEM image of a picked-up thin film

## Metal deposition unit

By installing the metal deposition unit, it is possible to deposit a carbon film on a desired position for protecting this position. The cartridge for a gas source is replaceable from the outside of a vacuum.



Metal deposition unit

# Specifications and Optional Attachments

## Principal specifications

Ion source	Ga liquid metal ion source
Accelerating voltage	5 to 30 kV (in 5 kV steps)
Magnification	×50 (for searching field) ×150 to ×300,000
Image resolution	6 nm (at 30 kV)
Maximum beam current	30 nA (at 30 kV)
Ion beam shapes during milling	Rectangle, line and spot
Specimen stage	Side-entry goniometer stage for TEM specimens
Vacuum pumps	SIP/TMP/RP

## Optional attachments

### ● Side Entry Goniometer Stage for TEM Specimens

EM-0261FSEG

### ● FIB Tip-on Holder

EM-02210

### ● FIB Bulk-Specimen Holder

EM-02220

### ● Bulk-Specimen Motor Stage

EM-02550FBSS

### ● FIB Bulk-Specimen Holder 1

EM-02560FBSh1

### ● FIB Bulk-Specimen Holder 2

EM-02570FBSh2

### ● FE-SEM Specimen Holder Adapter

EM-02580FSha

### ● Airlock System

EM-02590FALS

### ● Probe Current Detector

EM-02620FPCD

### ● Metal Deposition Unit

EM-02630FMDU

### ● Automatic Processing Software

EM-02520APS

### ● Specimen Pickup System

EM-02230

### ● Shuttle Retainer

EM-02280

### ● Lamp Holder

EM-02990

### ● Operation Keyboard


EM-02640FOKB

## Installation Requirements


Power	Single phase 100 V $\pm$ 10%, 50/60 Hz, 3 kVA
Grounding terminal	One, 100 $\Omega$ or less
Dry nitrogen gas	0.4 to 0.7 MPa
Exhaust duct	Exhaust pipe of 25 mm I. D.
Room temperature	18 to 25°C
Humidity	60% or less
Stray magnetic field	5 $\mu$ T or less
Floor vibration	0.5 $\mu$ m (p-p) or less
Noise	65 dB or less in the range up to 300 Hz

\*Dry nitrogen gas and a connecting tube should be provided by customer.

## Installation Room Example



\*Specifications subject to change without notice.



Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law" of Japan in compliance with international security export control. JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

**JEOL** JEOL Ltd. 1-2 Musashino 3-chome Akishima Tokyo 196-8558 Japan Sales Division ☎(042)528-3381 ㊚(042)528-3386 <http://www.jeol.co.jp/>

**ARGENTINA**  
COASIN S. A. C. I. yF.  
Virrey del Pino 4071,  
1430 Buenos Aires, Argentina  
Telephone: 54-11-4552-3185  
Facsimile: 54-11-4555-3321

**AUSTRALIA & NEW ZEALAND**  
JEOL (AUSTRALASIA) Pty. Ltd.  
Unit 9, 750-752 Pittwater Road,  
Brookvale, N. S. W. 2100, Australia  
Telephone: 61-2-9905-8255  
Facsimile: 61-2-9905-8286

**AUSTRIA**  
LABCO GmbH  
Dr.-Trittmel-Gasse 8,  
A-3013 Pressbaum, Austria  
Telephone: 43-2233-53838  
Facsimile: 43-2233-53176

**BANGLADESH**  
A.Q. CHOWDHURY & CO. Pvt. Ltd.  
Baridhara Central Plaza 87, Suhrawardy Avenue  
2nd Floor Baridhara, Dhaka-12129 Bangladesh  
Telephone: 880-2-9862272, 9894583  
Facsimile: 880-2-988070

**BELGIUM**  
JEOL (EUROPE) B. V.  
Planet II, Building B Leuvensesteenweg 542,  
B-1930 Zaventem, Belgium  
Telephone: 32-2-720-0560  
Facsimile: 32-2-720-6134

**BRAZIL**  
FUGIWARA ENTERPRISES  
INSTRUMENTOS CIENTIFICOS LTDA.  
Avenida Itaberaba,3563  
02739-000 Sao Paulo, SPI Brazil  
Telephone: 55-11-3983-8144  
Facsimile: 55-11-3983-8140

**CANADA**  
JEOL CANADA, INC.  
(Represented by Soquelec, Ltd.)  
5757 Cavendish Boulevard, Suite 540,  
Montreal, Quebec H4W 2W6, Canada  
Telephone: 1-514-482-6427  
Facsimile: 1-514-482-1929

**CHILE**  
TECSIS LTDA.  
Avenida Holanda 1248,  
Casilla 50/9 Correo 9, Providencia, Santiago, Chile  
Telephone: 56-2-205-1313  
Facsimile: 56-2-225-0759

**CHINA**  
JEOL LTD., BEIJING OFFICE  
Room No. B2308, Vantone New World Plaza, No. 2 Fuwai  
Street, Xicheng District, Beijing 100037, P. R. China  
Telephone: 86-10-6804-6321/6322/6323  
Facsimile: 86-10-6804-6324

JEOL LTD., SHANGHAI OFFICE  
Sanhe Building 11 F2, Yan Ping Road No. 121,  
Shanghai 200042, P. R. China  
Telephone: 86-21-6246-2353  
Facsimile: 86-21-6246-2636

JEOL LTD., GUANG ZHOU OFFICE  
S2204 World Trade Center Building #371-375, Huan Shi  
East-Road, Guangzhou, P. R. China  
Telephone: 86-20-8778-7848  
Facsimile: 86-20-8778-4268

JEOL LTD., CHENGDU OFFICE  
1807A Zongfu Bld., No. 45, Zhongfu Road Chengdu,  
Sichuan, P. R. China  
Telephone: 86-28-8662-2554  
Facsimile: 86-28-8662-2564

JEOL LTD., WUHAN OFFICE  
Room. 3216, World Trading Bld., 686 Jiefang Street,  
Hankou, Wuhan, P. R. China  
Telephone: 86-27-8544-8953  
Facsimile: 86-27-8544-8695

**FARMING LTD.**  
Unit 1009, 10/F., MLC Millennia Plaza  
663 King's Road, North Point, Hong Kong  
Telephone: 852-2815-7299  
Facsimile: 852-2581-4635

**CYPRUS**  
MESLO LTD.  
Scientific & Laboratory Division,  
P. O. Box 27709, Nicosia Cyprus  
Telephone: 357-2-666070  
Facsimile: 357-2-660355

**EGYPT**  
JEOL SERVICE BUREAU  
3rd Fl. Nile Center Bldg.,  
Nawal Street, Dokki, (Cairo), Egypt  
Telephone: 20-2-335-7220  
Facsimile: 20-2-338-4186

**FRANCE**  
JEOL (EUROPE) SAS  
Espace Claude Monet,  
1 Allée de Giverny 78290  
Croissy-sur-Seine, France  
Telephone: 33-13015-3737  
Facsimile: 33-13015-3747

**GERMANY**  
JEOL (GERMANY) GmbH  
Oskar-Von-Miller-Strasse 1,  
85386 Eching Germany  
Telephone: 49-8165-77346  
Facsimile: 49-8165-77512

**GREAT BRITAIN & IRELAND**  
JEOL (U.K.) LTD.  
JEOL House, Silver Court, Watchmead,  
Welwyn, Garden City, Herts AL7 1LT., England  
Telephone: 44-1707-377117  
Facsimile: 44-1707-373254

**GREECE**  
N. ASTERIDIS S. A.  
56-58, S. Trikoupi Str. P.O.Box 26140  
GR-10022 Athens, Greece  
Telephone: 30-1-823-5383  
Facsimile: 30-1-823-9567

**INDIA**  
Blue Star LTD. (HQ)  
Analytical Instruments Department,  
'Sahas'414/2 Veer Savarkar Marg,  
Prabhadary Mumbai 400 025, India  
Telephone: 91-22-5666-4068  
Facsimile: 91-22-5666-4001

Blue Star LTD. (New Delhi)  
Analytical Instruments Department,  
E-44/12 Okhla Industrial Area,  
Phase-11, New Delhi 110 020, India  
Telephone: 91-11-5149-4000  
Facsimile: 91-11-5149-4004

Blue Star LTD. (Calcutta)  
Analytical Instruments Department,  
7, Hare Street Calcutta 700 001, India  
Telephone: 91-33-2248-0131  
Facsimile: 91-33-2248-1599

Blue Star LTD. (Chennai)  
Analytical Instruments Department,  
Garuda Building 46, Cathedral Road  
Chennai 600 086, India  
Telephone: 91-44-5244-7210  
Facsimile: 91-44-5244-4190

**INDONESIA**  
PT. TEKNOLABindo PENTA PERKASA  
J1. Gading Bukit Raya,  
Komplek Gading Bukit Indah Blok I/11,  
Kelapa Gading Jakarta 14240, Indonesia  
Telephone: 62-21-45847057/58/59  
Facsimile: 62-21-45842729

**IRAN**  
IMACO LTD.  
No. 141 Estepoon Ave.,  
P. O. Box 13145-537, Tehran, Iran  
Telephone: 98-21-6402191/6404148  
Facsimile: 98-21-8978164

**ITALY**  
JEOL (ITALIA) S.p.A.  
Centro Direzionale Green Office Via Dei Tulipani,  
1, 20090 Pieve, Emanuele (MI), Italy  
Telephone: 39-2-9041431  
Facsimile: 39-2-90414353

**KOREA**  
JEOL KOREA LTD.  
Sunmin Bldg. 6th F1., 218-16, Nonhyun-Dong,  
Kangnam-Ku, Seoul, 135-010, Korea  
Telephone: 82-2-511-5501  
Facsimile: 82-2-511-2635

**KUWAIT**  
YUSUF I. AL-GHANIM & CO. (YIACO)  
P. O. Box 435, 13005 - Safat, Kuwait  
Telephone: 965-4832600/4814358  
Facsimile: 965-4844954/4833612

**MALAYSIA**  
JEOL (MALAYSIA) SDN. BHD. (359011-M)  
205, Block A, Mezzanine Floor,  
Kelana Business Center 97,  
Jalan SS 7/2, Kelana Jaya,  
47301 Petaling Jaya, Selangor, Malaysia  
Telephone: 60-3-7492-7722  
Facsimile: 60-3-7492-7723

**MEXICO**  
JEOL DE MEXICO S.A. DE C.V.  
Av. Amsterdam #46 DEPS. 402  
Col. Hipodromo, 06100 Mexico D.F. Mexico  
Telephone: 52-5-55-211-4511  
Facsimile: 52-5-55-211-0720

**PAKISTAN**  
Analytical Measuring System (Pvt.) Limited. AMS House  
Plot # 14C, Main Sahar Commercial Avenue,  
Commercial Lane 4,  
Khayaban-e-Sehar,  
D.H.A Phase 7, Karachi, Pakistan  
Telephone: 92-21-5345581/5340747  
Facsimile: 92-21-5345582

**PANAMA**  
PROMED S.A.  
Parque Industrial Costa del Este  
Urbanizacion Costa del Este  
Apartado 6281, Panama, Panama  
Telephone: 507-269-0044  
Facsimile: 507-263-5622

**PHILIPPINES**  
PHILAB INDUSTRIES INC.  
7487 Bagtikan Street, SAV Makati, 1203 Metro,  
Manila Philippines  
Telephone: 63-2-896-7218  
Facsimile: 63-2-897-7732

**PORTUGAL**  
Izasa. Portugal Lda.  
R. do Proletariado 1,  
2790-138 CARNAXIDE Portugal  
Telephone: 351-21-424-7300  
Facsimile: 351-21-418-6020

**SAUDI ARABIA**  
ABDULREHMAN ALGOSAIBI G. T.B.  
Algosai Bldg. Airport Rd., P. O. Box 215,  
Riyadh 11411, Saudi Arabia  
Telephone: 966-1-479-3000  
Facsimile: 966-1-477-1374

**SCANDINAVIA**  
JEOL (SKANDINAVISKA) A.B.  
Hammarbacken 6 A, Box 716  
191 27 Sollentuna, Sweden  
Telephone: 46-8-28-2800  
Facsimile: 46-8-29-1647

**SERVICE & INFORMATION OFFICE**  
JEOL NORWAY  
Ole Deviks vei 28, N-0614 Oslo, Norway  
Telephone: 47-2-2-64-7930  
Facsimile: 47-2-2-65-0619

**JEOL FINLAND**  
Ylakaupinkuja 2, FIN-02360 Espoo, Finland  
Telephone: 358-9-8129-0350  
Facsimile: 358-9-8129-0351

**JEOL DENMARK**  
Naverland 2, DK-2600 Glostrup, Denmark  
Telephone: 45-4345-3434  
Facsimile: 45-4345-3433

**SINGAPORE**  
JEOL ASIA PTE. LTD.  
29 International Business Park,  
#04-02A Acer Building,  
Tower B Singapore 609923  
Telephone: 65-6565-9989  
Facsimile: 65-6565-7552

**SOUTH AFRICA**  
ADI Scientific (Pty) Ltd.  
109 Blandford Road, North Riding, Randburg  
(PO box 71295 Bryanston 2021)  
Republic of South Africa  
Telephone: 27-11-462-1363  
Facsimile: 27-11-462-1466

**SPAIN**  
IZASA. S.A.  
Aragoneses, 13,  
28100 Alcobendas,  
(Poligono Industrial) Madrid, Spain  
Telephone: 34-91-663-0500  
Facsimile: 34-91-663-0545

**SWITZERLAND**  
JEOL (GERMANY) GmbH  
Oskar-Von-Miller Strasse 1,  
85386 Eching Germany  
Telephone: 49-8165-77346  
Facsimile: 49-8165-77512

**TAIWAN**  
JIE DONG CO., LTD.  
7F, 112, Chung Hsiao East Road, Section 1, Taipei,  
Taiwan 10023, Republic of China  
Telephone: 886-2-2395-2978  
Facsimile: 886-2-2322-4655

**JEOL TAIWAN SEMICONDUCTORS LTD.**  
11F, No. 346, Pei-Ta Road, Hsin-Chu City 300,  
Taiwan Republic of China  
Telephone: 886-3-523-8490  
Facsimile: 886-2-523-8503

**THAILAND**  
BEOTHAJ BANGKOK EQUIPMENT  
& CHEMICAL CO., Ltd.  
300 Phaholyothin Rd. Phayathai,  
Bangkok 10400, Thailand  
Telephone: 66-2-615-2929  
Facsimile: 66-2-615-2350/2351

**THE NETHERLANDS**  
JEOL (EUROPE) B.V.  
Tupolevslaan 28-A, 1119 NZ Schiphol-Rijk,  
The Netherlands  
Telephone: 31-20-6533088  
Facsimile: 31-20-6531328

**TURKEY**  
TEKSER LTD. STI.  
Acibadem Cad. Erdem Sok. Bayer Art. 6/1  
34660 Uskudar/Istanbul-Turkey  
Telephone: 90-216-3274041  
Facsimile: 90-216-3274046

**UAE**  
BUSINESS COMMUNICATIONS LLC.  
P. O. Box 2534, Abu Dhabi UAE  
Telephone: 971-2-6348495  
Facsimile: 971-2-6316465

**USA**  
JEOL USA, INC.  
11 Dearborn Road, Peabody, MA. 01960, U. S. A.  
Telephone: 1-978-535-5900  
Facsimile: 1-978-536-2205/2206

**JEQ USA, INC. WEST OFFICE**  
5653 Stoneridge Drive Suite  
#110 Pleasanton, CA. 94588 U. S. A.  
Tel: 1-925-737-1740  
Fax: 1-925-737-1749

**VENEZUELA**  
MITSUBISHI VENEZOLANA C. A.  
Avenida Francisco de Miranda Edificio Parque Canaima,  
Piso 2 Los Palos Grandes, Caracas, Venezuela  
Telephone: 58-212-209-7402  
Facsimile: 58-212-209-7496