

JEOL

GCmate, LCmate, RSVP

Chromatographic Enhancement Software

Component Detection for GC/MS and LC/MS

Introduction

A high chemical background or the continuous introduction of a reference compound can make it difficult to locate and identify trace components in a total ion current chromatogram in a GC/MS or LC/MS experiment. The *Shrader System for Windows™* software¹ available with the **JEOL GCmate, LCmate, and RSVP** mass spectrometers provides a solution to this problem based on the *Component Detection Algorithm (CODA)*, developed by Windig *et. al.*² The CODA algorithm identifies significant peaks by comparing a smoothed, mean-subtracted reconstructed ion chromatogram (RIC) with the original RIC.

Figure 1 shows the data processing screen when CODA is applied to a GC/MS chromatogram where a high chemical background is present. Note the improvement in signal-to-noise for the CODA-enhanced chromatogram. This enhancement can be dramatic and can often reveal "hidden" features in a chromatogram (see Figure 4).




Figure 1. CODA is used to enhance the signal-to-noise for a GC/MS reconstructed ion chromatogram (JEOL GCmate).




Figure 2. CODA is used to find series of related compounds in the mixture shown in Figure 1.

CODA can identify groups of mass-to-charge ratios whose mass chromatograms have similar elution times. This can be helpful in identifying related compounds. Figure 2 shows how the CODA-reconstructed chromatogram in Figure 1 can be decomposed into components.




Figure 3. Component detection for tea analyzed by electrospray ionization LC/MS (JEOL LCmate).

Figure 4 shows the result of applying CODA to a GC/MS exact-mass measurement where the reference compound (PFK) is continuously introduced, producing a high chemical background that makes it difficult to find the analyte peak.




Figure 4. CODA is used to eliminate background from continuous introduction of a reference compound (PFK) during an exact mass measurement with the JEOL GCmate.

References

¹ Shrader Analytical & Consulting Laboratories Inc., 3814 Vinewood, Detroit, MI 48208

² W. Windig, J.M. Phalp and A.W. Payne, **Analytical Chemistry**, 68(20), p. 3602.

JEOL USA, Inc.
 11 Dearborn Rd., Peabody, MA 01960
 Tel: (978)-535-5900 Fax: (978)-536-2205
 e-mail: ms@jeol.com WWW: <http://www.jeol.com>