



AccuTOF-GCx Series

Introduction of latest GC/HR-TOFMS system: JMS-T200GC AccuTOF-GCx - High Mass Accuracy and Stability

Introduction

JEOL has recently announced a 4th generation GC/HR-TOFMS system, the JMS-T200GC AccuTOF-GCx, in 2015. The AccuTOF-GCx offers high sensitivity, high mass resolving power, high mass accuracy, and a wide dynamic range in combination with high-speed data acquisition.

In this application note, we show we show high mass accuracy and m/z stability using this latest GC/HR-TOFMS system.

Results and Discussion

We made 100 replicate splitless injections of a 100fg/ μ L solution of octafluoronaphthalene (OFN) and acquired the mass spectra in electron ionization (EI) mode. All data were mass-calibrated with a one-point drift compensation by using a background column-bleed ion from the capillary column (m/z 281.05114: $C_7H_{21}O_4Si_4^+$) as a lock mass.

We obtained excellent reproducibility for the mass accuracy for the OFN molecular ion ($C_{10}F_8^+$, m/z 271.9867). All 100 measurements were within ± 1 mDa using the column-bleed drift compensation (Figure 2).




Figure 2. Mass error of the OFN 100fg/ μ L molecular ion ($C_{10}F_8^+$, m/z 271.9867) during 100 replicate measurements



Figure 1. 4th generation of the JEOL GC/HR-TOFMS system JMS-T200GC AccuTOF-GCx

Conclusion

The AccuTOF-GCx system provides high mass accuracy and stability by using a simple one-point drift compensation.




Figure 3. Overlay chart of the OFN molecular ion peaks