

Scientific / Metrology Instruments CROSS SECTION POLISHER[™]

Solutions for Innovation

CROSS SECTION POLISHER[™]

IB-10500HMS High Throughput Milling System





High milling rate of cross-section achieved by the new ion source: **1.2 mm/h or more***² (**2.4 times** than the previous milling rate.)

The high throughput milling system optimizes the ion source electrodes and enables higher accelerating voltages, thus improving the ion-beam current density. Our newly developed ion source achieves a high milling rate of cross-section of 1.2 mm/h or more (2.4 times than the previous milling rate.)

Cross-section milling rate of the new ion source

Specimen: Silicon wafer, Accelerating voltage: 10 kV, Milling time: 1 h





Cross-section milling of powders



Accelerating voltage: 10 kV, Milling time: 10 min

The specimen can be processed in a short time. This leads to enhancement of work efficiency.

Cross-section milling of a low melting-point alloy (cooling)

Accelerating voltage: 10 kV, Milling time: 30 min

The right SEM image shows an Sn-Bi alloy with a melting point of 150 °C. A low melting-point metal can be melted due to the processing heat; therefore, cooling of the metal is required before milling. High throughput milling is applied to the heat-sensitive specimen while the specimen is kept cooled *3 . Then, a cross-section specimen with a reduced heat damage is obtained in a short time.





10 µm



Large Area Milling*1, *4

Planar Surface Milling of Larger Area

The new high throughput milling system has enabled the irradiation of an ion beam onto a larger area of the specimen. Planar surface milling is effective to remove scratches generated on the specimen surface or crystalline strains, which are caused by mechanical polishing.

Planar surface milling of a concrete

Accelerating voltage: 10 kV, Milling time: 20 min

Large-area planar surface milling was applied to a concrete with a width of 20 mm. After milling, polish scratches and contamination were removed, allowing for clear observation of particles of stone and cement contained in a concrete.



*1 This function is included in IB-19530CP or IB-19520CCP which incorporates the optional high throughput specification. *2 Milling of 1 h, Si equivalent, Edge distance: 100 µm *3 This function is included in IB-19520CCP. *4 Large Specimen Rotation Holder, IB-11550LSRH, is required.

Specifications		Standard specification		High throughput specification *1 *2		
		IB-19530CP	IB-19520CCP	IB-19530CP + IB-10500HMS	IB-19520CCP + IB-10500HMS	
Ion accelerating voltage		2 to 8 kV		2 to 10 kV		
Milling speed		500 $\mu m/h$ or more (accelerating voltage 8 kV) \star_3		1200 $\mu m/h$ or more (accelerating voltage 10 kV) $^{\star 4}$		
Specimen swing function *5		Auto specimen swing by $\pm 30^{\circ}$		Auto specimen swing by ± 30°, Angle setting swing		
Auto milling start mode		0	0	0	0	
Auto cooling milling start mode / Auto return to room temperature mode		-	0	-	0	
Specimen stage ultimate cooling temperature		-	–120 °C or less	-	–120 °C or less	
Cooling temperature settable range		-	–120 to 0 °C	-	–120 to 0 °C	
Specimen cooling time to reach –100 °C		-	Within 60 min	-	Within 60 min	
Specimen cooling retention time		-	8 h or more *6	-	8 h or more *6	
Air isolation function		-	0	-	-	
Intermittent milling mode		Ion beam irradiation time and stop time are settable (ON: 1 to 999 s, OFF: 1 to 999 s)				
Fine milling mode		Milling conditions automatically switched				
Large-area cross-section milling mode *7		Maximum milling width: 8 mm (with optional Large Area Milling Holder IB-11730LMH)				
Large-area planar surface milling mode		-	-	O *8	○ *8	
Maximum specimen size	Cross-section milling	11 mm (W) × 10 mm (L) × 2 mm (T) (with standard holder for IB-19530CP) 11 mm (W) × 8 mm (L) × 3 mm (T) (with standard holder for IB-19520CCP) 25 mm (W) × 15 mm (L) × 10 mm (T) (with optional Large Area Milling Holder IB-11730LMH)				
	Planar surface milling	40 mm (diameter) × 15 mm (T) (with optional Large Specimen Rotation Holder IB-11550LSRH)				
Specimen movements		X-axis: ± 6 mm, Y-axis: ± 2.5 mm				
Operation		Touch panel, 6.5-inch display				
Positioning for milling		Monitor from above the specimen stage with a camera *9. Milling position is adjustable with an optical microscope.				
Positioning camera (magnification)		Approx. ×70 (on 6.5-inch display)				
Monitoring camera (magnification)		Approx. ×20 to 100 (on 6.5-inch display) *Note: When used with IB-19530CP + IB-14510MCAM *10 or IB-19520CCP.				
External monitor output		Positioning camera and Monitoring camera can be switched for displaying one on the external monitor. *Note: When used with IB-19530CP + IB-14510MCAM * ¹⁰ or IB-19520CCP + EC-10020VST * ¹¹ .				
Preset function		4 sets of milling conditions (accelerating voltage, Ar gas flow, milling time, intermittent milling)				
Dimensions and weights: Basic unit		545 mm (W) × 550 mm (D) × 420 mm (H), Approx. 66 kg (with IB-19530CP + IB-14510MCAM attached) 690 mm (W) × 720 mm (D) × 530 mm (H), Approx. 75 kg (with IB-19520CCP attached)				
Dimensions and w	eights: Rotary pump	1	50 mm (W) × 427 mm (D) :	× 230 mm (H), Approx. 16 k	g	

Installation Requirements

Power supply	Single phase 100 to 120 V AC, 50/60 Hz, Allowable input voltage fluctuation: less than 10%, Rating: 15 A or more			
Maximum power consumption	650 VA			
Grounding	100 Ω or less			
Argon gas *12	Dry argon, Purity: 99.9999% or more. Pressure: 0.1 to 0.2 MPa (1.0 to 2.0 kgf/cm ²), Hose joint: ISO 7/1 Rc 1/4			
Room temperature	15 to 25 °C			
Room humidity	60% or less (no condensation)			

This is optional, which is added at the time of shipment from factory.

*2 An ion source and an ion-current detection unit are different for the standard specification and the high throughput specification.

Average over 2 h, Si equivalent, Edge distance 100 μm Milling of 1 h, Si equivalent, Edge distance 100 μm Patent No. (Japan): 4557130 *3

*4

*5

- *6
- As the set temperature is higher, the cooling retention time is longer. This mode can be used in combination with the cooling function of IB-19520CCP.

*8 When used with IB-11550LSRH. *9 Patent No. (Japan): 4208658

JEOL

*10 With IB-14510MCAM attached, the specimen can be monitored in real time.

JEOL Ltd.

- The status of the specimen can be observed while milling is in progress.
- The external monitor must be prepared by the customer. *11 With EC-10020VST attached, the camera image can be displayed on the external monitor. The external monitor must be prepared by the customer.
- *The specifications and appearance of the instrument are subject to change without notice. Certain products in this brochure are controlled under the "Foreign Exchange and Foreign Trade Law' of Japan in compliance with international security export control, JEOL Ltd. must provide the Japanese Government with "End-user's Statement of Assurance" and "End-use Certificate" in order to obtain the export license needed for export from Japan. If the product to be exported is in this category, the end user will be asked to fill in these certificate forms.

*The screen images in the catalog include items that are still under development,

*12 The argon gas, gas cylinders and regulator must be prepared by the customer.



AUSTRALIA & NEW ZEALAND /JEOL (AUSTRALASIA) Pty.Ltd. Suite 1, L2 18 Aquatic Drive - Frenchs Forest NSW 2086 Australia • BELGIUM /JEOL (EUROPE) B.V. Planet II, Gebouw B Leuvensesteenweg 542, B-1930 Zaventem Belgium
BRAZIL /JEOL Brasil Instrumentos Cientificos Ltda, Av, Jabaquara, 2958 5° andar conjunto 52; 04046-500 Sao Paulo, SP Brazil • CANADA /JEOL CANADA, INC. 3275 tere Rue, Local #8 SHubert, QC J3Y-8Y6, Canada • CHINA /JEOL (BELJING)
Co., LTD. Zhongkezjuvan Building South Tower 2F, Zhongguancun Nansanjis Street No. 6, Haidan District, Beijing, P.R.China • EGYPT /JEOL SERVICE BUREAU 3rd FI. Nile Center Bldg., Nawal Street, Dokki, (Cario), Egypt • FRANCE /JEOL
(EUROPE) SAS Espace Claude Monet, 1 Allee de Giverny 72820, Croisey-suc-Seine, France • GERMANY) GmbH Gute Aenger 30 85366 Freising, Germany • GREAT BRITIN & IRELAND /JEOL (UK) LTJ. JEOL House, Silver
Court, Watchmead, Welwyn Garden City, Herts AL7 1LT, U.K. • INDIA /JEOL INDIA PVT, LTD. Unit No.305, 3rd Floor, ABW Elegance Tower, Jasola District Centre, New Delhi 110 025, India /JEOL INDIA PVT, LTD. Hyderabad office 422, Regus
Soliatire Business centre. 1-10-39 to 44, Jevel 4, Gumidelli Towers, Old Airport Road, Begumpet Hyderabad - 500016, India • ITALY /JEOL (ITALIA) S.p.A. Palazzo Pacinotti - Milano 3 City, Via Ludovico IMoro, 6/A 2008 Basiglio(MI) Ital) • KOREA
JEOL KOREA LTD. Dongwoo Bidg, TF, 1443, Yangia Baero, Gangdong-Gu, Seou), G5355, Korea • MALAYSIA /JEOL (UK) LCMLAVSIA) SDN.BHD. 508, Block A, Level 5, Kelana Business Centre, 97, Jalan SS 7/2, Kelana Jaya, 47301 Petaling Jaya,
Selangor, Malaysia • MEXICO /JEOL DE MEXICO S.A. DE C.V. Arkansas 11 Piso 2 Colonia Napoles Delegacion Benito Juarez, C.P. 03810 Mexico D.F., Mexico • QATAR Mannai Trading Company V.L.L. ALI Emadi. Complex, Salwa Road P.O. Box
Pr. Doroporation Road #01-12 Corporation Road #01-12 Corporation Road #04-012 Corporation Road #04-05192 Science 4,
NL2153 PH Nieuw-Vennep, The Netherlands • USA /JEOL USA, INC 11 Dearborn Road, Peabody,

and are subject to change without notice.