

Key Word: AES, Attachmnets, Option Holders

# Sample Holders for JAMP series

The JAMP series (Auger Microprobe) has 7 types of holder for various samples. These are prepared for any sizes, any shape, and any analysis conditions. Please contact to JEOL Ltd. if you need a custom-made holder !

<p>Standard 12 mm holder</p>  <p>Parts No. 8043 97074</p>	<p>Standard 20 mm holder</p>  <p>Parts No. 8043 97066</p>	<p>25 mm holder (1 inch size)</p>  <p>Parts No. 7801 63818</p>
--	--	---

Sample Size: 12 mm $\Phi$  × 4 mm

Sample Size: 20 mm $\Phi$  × 5 mm

Sample Size: 26 mm $\Phi$  × 7 mm

Attention: The sample is not a JEOL products in right side

<p>Quick Holder</p>  <p>Parts No. 8044 13223</p>	<p>45 deg. Pre-tilt Holder (for sample cut by FIB)</p>  <p>Parts No. 8044 08670</p>	<p>50 mm Holder</p>  <p>Parts No. 8044 04577</p>
---	--	---


Sample Size: 12 mm $\Phi$  × 4 mm

Sample Size: 12 mm $\Phi$  × 4 mm

Sample Size: 50 mm $\Phi$  × 6 mm

Attention: the standard stage unit can allow an analysis area within 20 mm × 20 mm in the holder center.


Sample Holder for irregular shape



Parts No. 8044 05964

Sample Size: 26 mm × 10 mm × 7 mm(t)

Sample Holder for cross-section sample



Parts No. 8044 06367

Sample Size: 12 mm × 12 mm × 7 mm (t)

Welcome your order of a custom-made holder if you need. Please contact us !

Access the QR codes below for more information on the Field Emission Auger Microprobe

◆ Overview →



◆ Mechanisms →

