



Perfluorotributylamine (PFTBA) (EI, positive-ion)

Disclaimer: These data are for informational use only. No guarantees are made for the accuracy of the data presented.

Introduced via: Batch reference inlet

Exact Mass	Composition
68.99521	CF ₃
99.993613	C ₂ F ₄
118.992016	C ₂ F ₅
130.992016	C ₃ F ₅
218.985629	C ₄ F ₉
263.9871	C ₅ F ₁₀ N
313.983917	C ₆ F ₁₂ N
325.983917	C ₇ F ₁₂ N
375.980743	C ₈ F ₁₄ N
413.977539	C ₈ F ₁₆ N
425.977539	C ₉ F ₁₆ N
463.974335	C ₉ F ₁₈ N
501.971161	C ₉ F ₂₀ N
575.967957	C ₁₂ F ₂₂ N
613.964722	C ₁₂ F ₂₄ N